911 NETFULL

880.02

CE
Regulation (EU)
2017/745

EN

WWW.KONG.IT

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INDEX

	1 - SYMBOLS AND ASSISTANCE	4
	1.1 Symbols	4
	1.2 Assistance	4
	2 - GENERAL INFORMATION	5
	3 - TECHNICAL CHARACTERISTICS	6
	3.1 Terminology and materials of components	6
	3.2 Dimensions	7
	3.3 Capacity	7
	3.4 Optional components and spare parts	8
	4 - SPECIFIC INFORMATION	10
	4.1 Intended use	10
	4.2 Assembly	10
	4.3 Positioning and physically securing patient	12
	4.4 Transport of patient	13
	5 - MAINTENANCE AND REPAIR	15
	5.1 General	15
	5.2 Maintenance	15
	5.3 Repair	15
	6 - STORAGE	16
	7 - CHECKS, INSPECTIONS AND OVERHAULS	17
	7.1 Checks	17
	7.2 Inspections	17
	7.3 Overhauls	17
	8 - DURATION AND WARRANTY	18
	8.1 Product Lifetime	18
	8.2 Disposal	18
	8.3 Warranty	18
	8.4 Legal obligations	18
	9 - LABELING AND SYMBOLS	19
	9.1 Medical Device label	19
	9.2 Symbols	19
	10 - DOCUMENTS	20
	10.1 Inspections and servicing register	20
	10.2 Maintenance and repair register	21
	10.3 EU Declaration of conformity	22
	11 - LIST OF STANDARDS AND REFERENCE STANDARDS APPLIED	23
	11.1 Standards applied	23
	11.2 Standards used as a reference	24

1.1 SYMBOLS

To make the manual easy and clear to read, the following symbols are used for important warnings for correct, safe use of the device.

REQUIREMENT FOR CORRECT USE

Identifies the presence of information for correct use of the device.

INFORMATION REQUIREMENT

Identifies the presence of useful, general information, the reading of which guides the user towards using the device with care and / or performing specific actions.

It identifies that the Medical Device is manufactured, designed and produced in compliance with the provisions of the General Safety and Performance Requirements of Regulation (EU) 2017/745 (Class I medical device, in compliance with classification rule 1 as indicated in the attachment VIII).

1.2 ASSISTANCE

For information, contact the Kong Customer Service Department by:

telephone +39 0341 630506

- fax +39 0341 641550

- email: safetycare@kong.it

To facilitate the assistance we can give you, always communicate or indicate the serial number (SN) shown on the label affixed to the Medical Device.

CHAPTER 2

2 GENERAL INFORMATION

The information provided by the manufacturer (hereinafter referred to as information) must be read and clearly understood by the user prior to using the Medical Device (MD hereafter). The information regards the description of the features, performance, assembly, disassembly, maintenance, preservation, disinfection, etc. of the device. Even though the information offers tips on use, this information shall not be considered as a user manual under actual conditions of use.

WARNINGS AND LIMITATIONS OF USE:

- this device shall be strictly used by people who are physically fit, trained (instructed and taught) to use the device and with specific experience regarding moving the patient or, during training activities, by people under direct supervision of the trainers/ supervisors who guarantee the safety thereof,
- do not use the device before fully reading and understanding this user manual,
- strictly follow the manufacturer's information, improper use of the device is hazardous,
- modifying and/or repairing the device is strictly forbidden,
- all checks described in chapter 7 must be carried out prior to and after using the device. In case of any doubt on the efficiency of the device, the user must replace it immediately,
- non-compliant use, deformations, falls, wear, chemical contamination, exposing textile/plastic components/devices to temperature below -30°C or above

- +50°C and metal components/devices to temperatures exceeding 100°C are some examples of causes that can reduce, limit and end the life of the device,
- Prior to any rescue operation, be keen not to exceed the capacity indicated in paragraph 3.3,
- in order to reduce risks of exposure to / transmission of infectious diseases, clean and disinfect the device as indicated in chapter 5,
- improper use of the patient immobilisation systems can jeopardise the safety of the patient,
- always check the compatibility of the devices used alongside the device by consulting the relative manufacturer's information,
- use of spare parts or optional components different from the ones indicated in paragraph 3.4 can be hazardous,
- do not expose the device to sources of heat and at contact with chemical substances. Reduce exposure to direct sunlight as much as possible. At low temperatures and in the presence of humidity, formation of ice could reduce flexibility and increase cutting and abrasion-related risks on textile and synthetic devices,
- Report to the manufacturer and the competent authority of the member state where the user and / or patient is established any serious incident occurring in relation to the MD.

All our devices are tested/ checked piece by piece in compliance with the procedures laid down by the Quality System certified in accordance with the UNI EN ISO 9001 standard. Laboratory tests, testing, information and standards do not always reproduce the practical result. Thus, the results obtained under the actual conditions of use of the device in the natural environment may differ, even considerably at times. The best information lies in the continuous practical use under the supervision of skilled/ expert/qualified people.

Warning: not suitable for use in an ATEX environment (Directive 94/9/CE)

CHAPTER 3

TECHNICAL FEATURES

3.1 PARTS TERMINOLOGY AND MATERIALS

- A – Aluminum alloy frame, removable in two symmetrical parts,
- B – Polypropylene net,
- C – Polyurethane foam mats,
- D – Polyester webbing for physically securing patient,
- E – Polyester tapes for fixing the mats,
- F – Aluminum alloy handles,
- G – Polyurethane padding for handles.


Fig.1

3.2 DIMENSIONS

Length without handles:	220 cm
Width with handles:	350 cm
Width:	60 cm
Total weight of parts (fig.1):	8 kg

3.3 CAPACITY

The table below defines the capacities of the “911 NETFULL” stretcher, based on the method of use deemed most appropriate by the rescue operator, the outcome of the laboratory tests and the safety coefficients.

How to use	Capacity
Transport with the handles (fig. 2)	150 kg
Lifting and lowering with rescue devices (fig. 3) <i>with 880020100KK PHOENIX hanging kit (optional components)</i>	
Lifting and lowering using helicopters (fig. 3) <i>with 880020100KK PHOENIX hanging kit (optional components)</i>	100 kg
Overall buoyancy thrust (fig. 3) <i>with the kit 880020300KK FLOK-IT (L) fully submerged</i>	400 N


Fig.2

Before any recovery operation, make sure that the weight does not exceed the capacity defined above!

3.4 OPTIONAL COMPONENTS AND SPARE PARTS

The methods of connection and use of the optional components are defined in the relative user instructions.

3.4.1 Optional components

880020100KK PHOENIX

Suspension with coloured ends (red, head side - blue, feet side) and connectors with screw-equipped ring nut, suitable for lifting and lowering using rescue devices and helicopters.



Fig.3

880020300KK FLOK-IT

Floating kit for stretchers, made of three cylinders of different lengths (1 m for the sides - 0,4 m head side), each of one equipped with two fastening straps with quick release buckles.


Fig.4

880010300KK 911 BAG
Winch carriage bag.


Fig.5

3.4.2 Spare parts

- C – Mattresses,
- G – Padding for handles.


Fig.6

CHAPTER 4

INFORMATION SPECIFICATIONS

4.1 INTENDED USE


The “**911 NETFULL**” stretcher is a MD specially designed for the recovery and transport of patients, even if immobilized on “**X-TRIM**” spinal boards or on the “**VACUUM**” vacuum mat. All decisions regarding movement and immobilization of the patient, as well as the duration of these operations, the methods to be adopted and use in conjunction with other devices, must be taken and implemented by expert, trained personnel only.

The “**PHOENIX**” (I) hanging kit, not supplied, makes the stretcher capable of being lifted by helicopter.

4.2 ASSEMBLY

Assemble the two parts of the structure in the following steps:

- Rotate the two safety arms (fig. 7),
- Lift and rotate the two locking arms (fig. 8),

c) Connect the two parts of the structure using the pins on the frame (fig. 9),


Fig.9

d) Reposition the two locking arms in their respective brackets (fig. 10),

e) Reposition the two safety arms (fig.11) and check the tightness of the coupling,

f) Fully tighten the two ring nuts (Fig. 12).


Fig.10

Fig.11

Fig.12

4.3 POSITIONING AND SECURING THE PATIENT

With the stretcher positioned on a flat, stable surface:

- Unhook the buckles and open the securing straps (Fig. 13),
- Lay the patient in the stretcher, hook on the buckles and tension the straps until they secure the patient (Fig. 14).

This procedure must also be used for patients immobilized on the “X-TRIM” spinal board or on the “VACUUM” vacuum mat.

Avoid prolonged contact of the patient’s skin with the fabrics of the device


Fig.13


Fig.14

4.4 PATIENT TRANSPORT

The stretcher enables the patient to be lifted and transported by:

- a) the side slots in the frame
- b) the handles,
- c) the suspension assembly.

4.4.1 Connecting the handles

- a) Connect the handles to the frame of the stretcher (fig.15), making sure to observe the coupling numbering,


Fig.15

- b) Lock the connection by inserting the shaped plug into the corresponding hole (Fig. 16).


Fig.16

4.4.2 Connection to the hanging kit

- a) Connect the hanging kit to the corresponding side holes on the stretcher, using the eight connectors (red, head side - blue, feet side) (Fig. 17),
- b) Screw the ring nuts of the connectors fully tight.


Fig.17

CHAPTER 5

MAINTENANCE AND REPAIR

5.1 GENERAL

The “**NETFULL**” stretcher is made from materials with high resistance to wear and to external agents. Despite this, practical conditions of use will require some maintenance to be carried out and, in certain special cases, even some repairs.

Maintenance and repairs must be recorded in an appropriate log, an example of which is given in chapter 10.

5.2 MAINTENANCE

The maintenance procedures that must be performed by the user are:

- a) cleaning: wash after each use with lukewarm drinking water (max 40 ° C), possibly with the addition of a neutral detergent (e.g. Marseille soap). Rinse and leave to dry in the shade, away from direct heat,
- b) disinfection, when deemed necessary or after prolonged inactivity (more than 21 days), preceded by cleaning: wipe with a clean cloth , soaked in a water solution of sodium hypochlorite (bleach) with a 0,1% concentration (1000 ppm). When blood or other bodily fluids are present, the recommended concentration of sodium hypochlorite is 0,5%(5000ppm).

5.3 REPAIR

Repairs may only be carried out by the manufacturer.
The user is permitted to replace with new, original parts only those parts mentioned in paragraph 3.4.2.

CHAPTER 6

STORAGE

6 STORAGE

After cleaning, disinfection and drying, store the device and its accessories in a dry (relative humidity 40-90%), cool (temperature 5-40 ° C), dark (avoid UV radiation), chemically neutral (absolutely avoid saline and / or acidic environments) place, away from sharp edges, sources of heat, humidity, corrosive substances or other possible damaging conditions.

Do not store this device wet!

7.1 CONTROLS

Before and after each use, it is necessary to check the device and to make sure that:

- it is suitable for the use for which it is designed,
- it has not undergone any mechanical deformations and shows no cracks or signs of wear,
- the textile parts do not show cuts, burns, residues of chemical products, excessive hair, wear; in particular, check the areas in contact with the metal components,
- all markings, including the labels, are legible.

7.2 INSPECTIONS

The device must be inspected at least once a year, starting from the date of first use, by competent personnel trained and authorized by KONG SpA.

The date of first use and the results of the inspections must be recorded in the Inspection and Overhauls Log, an example of which is given in chapter 10.

7.3 OVERHAULS

All overhauls of the device must be performed by KONG SpA, or by competent personnel trained and authorized by the same, if:

- any malfunctions are detected,
- the outcome of pre and post-use checks or inspections is negative.

The results of the overhauls must be recorded in the Inspections and Overhauls Log, an example of which is given in chapter 10.

The overhauled device is guaranteed for one year from the date of the overhaul

CHAPTER 8

PRODUCT LIFETIME AND WARRANTY

8.1 PRODUCT LIFETIME

Carefully read the section “Warnings and Limitations of Use” in chapter 2.

The lifetime of this device is determined by the results of inspections, carried out at least once a year, (paragraph 7.2) and of any overhauls.

The suspension assembly, even if never used, has a maximum lifetime of 10 years from the date of manufacture shown on the label.

Important: remove and make devices unusable that have not passed the pre-use, post-use or periodic inspections.

8.2 DISPOSAL

For correct disposal, you must comply with the regulations in force in the country of use or the waste disposal procedures used by the reference hospital.

8.3 WARRANTY

The manufacturer guarantees the conformity of the device with the standards in force at the time of production. The warranty for defects is limited to defects in raw materials and manufacturing, does not include normal wear and tear, oxidation, damage caused by non-compliant use and / or competitions, incorrect maintenance, transport, preservation or storage, etc. The warranty will become void immediately if any changes to the device are made or if it is tampered with. Its validity corresponds to the legal warranty for the country in which the device was sold, from the date of sale by the manufacturer. Once this period has expired, no claim can be made against the manufacturer. Any request for repair or replacement under warranty must be accompanied by proof of purchase. If the defect is acknowledged, the manufacturer undertakes to repair or, at its discretion, replace the device or else refund the purchase price. In no case shall the manufacturer's liability extend to an amount higher than the purchase price of the device.


8.4 LEGAL OBLIGATIONS

Professional and leisure activities are often regulated by specific national laws that may impose limits on and / or obligations regarding the use of these devices. It is the user's duty to know and to apply such laws, as these may provide for restrictions different from those stated in this information notice.

CHAPTER 9

LABELLING AND SYMBOLS


9.1 MEDICAL DEVICE LABEL

Any changes to the positioning of the symbols does not change the content thereof.

9.2 SYMBOLS

Identification and traceability

	Medical Device
	Product identification code
	Unique Device Identifier
	Manufacturer's identification
	Unique serial number
	Country, year and month of production.

Data matrix UDI
 Readable with the Apps:
 iGepir - BarValid - Barcode Scanner

Information and safety

	Consult the user manual
	Compliant with Regulation (EU) 2017/745

10.3 EU DECLARATION OF CONFORMITY

The EU Declaration of Conformity is available for download at www.kong.it

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11.1 APPLIED STANDARDS

- EN 597-1: 201
Furniture - Assessment of the ignitability of mattresses and upholstered bed bases -Part 1: Ignition source smouldering cigarette
- EN 1865-1: 2010 + A1: 201
Patient handling equipment used in road ambulances - Part 1: General stretcher systems and patient handling equipment
- EN 62366-1: 2015 + A1: 202
Medical devices - Application of usability engineering to medical devices.
- EN ISO 846:201
Plastics - Evaluation of the action of microorganisms
- EN ISO 10993-1: 202
Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- EN ISO 10993-5: 200
Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
- EN ISO 10993-23: 202
Biological evaluation of medical devices - Part 23: Tests for irritation
- EN ISO 14971: 201
Medical devices - Application of risk management to medical devices.
- EN ISO 15223-1: 202
Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
- EN ISO 20417: 202
Medical devices - Information to be supplied by the manufacturer

11.2 STANDARDS USED AS REFERENCE

- EN ISO 13485: 2011
Medical devices - Quality management systems - Requirements for regulatory purposes
- EN 13718-1: 2014 + A1: 2021
Medical vehicles and their equipment - Air ambulances - Part 1: Requirements for medical devices used in air ambulances

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