

STRETCHERS ACCESSORIES

VACUUM 3

NEW

VACUUM MATTRESS | FOR INJURED IMMOBILIZATION

Vacuum mattress for wounded immobilization. Its dimensions fit LECCO, 911, Rolly and Kit Everest stretchers, keeping the immobilization of the injured also in vertical position.

There is no interference with X-ray or MRI it is a clear image, the only visible shadow is from the valve, but very distinct that is a round object, the same as if the patients have a shirt with buttons on. The new Polyester beads stays very clear, even after years of using the splints, not as the old Styrofoam beads collapse (mature) and change of the density fooling the X-ray technician with an image of a fracture.

The temperature ranging from +150 C to -70 C without any problem with the Polyurethane fabric.

It is recommended a maximum load of 200 kg.

The fabrics of the carrying cases is 1000 Denier Michigan Cordura Nylon with 10 mm zipper and reflective tape. The handles and straps are 35 mm Polyester webbing with 35 mm double adjusting multi colored buckles.

Double laminated Polyurethane with a scrim of nylon inside to make it have a maximum penetration, abrasion and rip strength. The bottom material is a heavy 400 denier (0.45 mm) and the inner material is also a double laminated Polyurethane with a weight of 200 D. (0.35 mm).

Polyurethane materials are very tough and resist Benzene, Paint thinner, Acetone and gasoline, even Chlorine, but if will fade the colors.

All products are tested for a 24 hrs vacuum hold before leaving the factory. It is a life time warranty on all welding seams and the valve plug. One year warranty on the handles and straps. CE Certification and FDA approval in USA.

Repair kit with Polyurethane glue and patches.

Vacuum pump and carriage bag included.



Recommended for

Professional: Rescue activities

Mountain rescue, Firebrigades and civil protection, Helirescue

Articolo	Ref.	DM	N. Reg. DM	CE EN	Material	▲ (g)	Colour
872050000KK	MT80111	Classe I	GMDN 35358	Direttiva 93/42/CEE	Polyamide	3500	Black Red